	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNING

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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The Duoning medium additive is a product prepared with dry powder Transpro CD 01 medium to prepare liquid Transpro CD 01. Transpro CD 01 is a universal transient medium that can be used for both HEK293 cells and CHO cells for passaging, high-density culture and transient transfection culture without the need for centrifugal exchange during transient transfection. Transpro CD 01 is suitable for transient transfection and expression culture of antibodies, recombinant proteins, and viruses during the research and development process using HEK293 series cells such as HEK 293, Expi293F, 293F, 293E, and CHO series cells such as ExpiCHOS and CHOS. This product is a fully chemically defined medium, no animal-derived components, no protein components, no animal or plant-derived protein hydrolysates, and no growth factors. This product does not contain HT or anti-clumping agents. Transpro CD 01 liquid package does not contain L-glutamine and requires additional supplementation of 4-6 mM L-glutamine for use. Transpro CD 01 medium powder package contains 6 mM L-glutamine.

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Taking Duoning medium additive combined with dry powder Transpro CD 01 medium to prepare Transpro CD 01 liquid medium as an example (1L for example)

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Prepare ultrapure water with a volume of about 90% (20~30°C);

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Add 0.165g of Duoning medium additive powder, stir for 5min;

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Stir for 30min until completely dissolved;

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Sterilized with 0.22 μ m filter.

【细胞培养 Cell culture】

① 建议细胞接种密度 Suggested cell inoculation density: 0.2~1.0 \times 10⁶ cells/mL

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Most cell lines use this product without any domestication, and can be directly inoculated into this medium and passed for more than three times. For some cell lines, domestication may be used when using this series of medium.

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Prepare frozen solution on the super clean workbench : 90% Transpro CD 01 +10% dimethyl sulfoxide (DMSO) mixed solution, precooling at 2~8°C (heat will be released when DMSO is diluted);

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Frozen cell fluid: seed cells were in the exponential growth period, the density is greater than 1.5×10^6 cells/mL, and the viability is greater than 95%.

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Cell fluid was centrifuged at 800rpm for 5 min;

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Slowly pour out the supernatant, resuspend the cells with cryopreservation solution, the cryopreservation density is $1.0 \sim 1.5 \times 10^7$ cells/mL, and transfer the cells to a sterile cryopreservation tube;

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Place the cryopreservation tube in the cryopreservation box containing isopropyl alcohol, freeze it at -80°C overnight, and then transfer it to the liquid nitrogen tank for long-term storage. If there is no freezing box, the temperature can be reduced manually by gradient as follows:

- 4°C 冻存 30min;
- freeze at 4°C for 30min;
- -20°C 冻存 2~4 小时;
- freeze at -20°C for 2~4h;
- -80°C 冻存过夜;
- freeze at -80°C overnight;
- 转移至液氮罐中长期贮存。
- transfer frozen cells to liquid nitrogen tank for long-term storage.

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Prepare 36.5°C warm water to thaw frozen cells;

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Take out the frozen tube from the liquid nitrogen tank and quickly thaw frozen cells in 36.5°C warm water;

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After wiping the cryopreservation tube with 75% ethanol, open the cryopreservation tube in the sterile operation table, transfer the cell fluid to a 15 ml centrifuge tube containing 2-5 mL of Transpro CD 01, blow and mix well, centrifuge at 800 rpm for 5 minutes;

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Slowly pour out the supernatant, resuspend with 15~20 ml preheated Transpro CD 01, and transfer to a 125 ml shake flask;

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Place it in a shaking incubator with 8% CO₂, 110 ~ 130rpm, at 36.5°C for culture;

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After 2~3 days of culture, the cells were counted and subcultured.

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The cells are seeded at 0.2~1.0×10⁶ cells/mL, count and subculture every 2 ~ 3 days. In the first three passages, the volume remained unchanged to restore cell viability. After the cell viability recovers to normal and reaches more than 90%. The seed cells were expanded at the density of 0.2~1.0×10⁶ cells/mL until reaching the required volume. The criteria for normal seed state: the viability was greater than 95%, the cell morphology was regular and round, and the growth doubling time was normal.

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Shanghai production base, powder packaging: 2°C to 8°C, protect from light; validity period: 24 months.

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Wuxi production base, powder packaging: 2°C to 8°C, protect from light; retest date: 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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